



LAMPIRAN


Lampiran 1 Binding pocket protein target 6XR3 dengan ligan uji




Kaempferol




Kuersetin



Lupeol




Luteolin




Warna Merah
Ligand


Lampiran 2 Binding pocket protein target 7JN2 dengan ligan uji




Asam Kafeat




Kaempferol



Kuersetin




Luteolin




Warna Merah
Ligand


Lampiran 3 Binding pocket protein target 7AAP dengan ligan uji




Asam Ferulat




Asam Kafeat



Kuersetin



Luteolin



Warna Merah
Ligand Asli